

## Physical Specifications

Dimensions (LXHXD)	177 X 252 X 202 mm
Weight	2.5 kg
Display	1.77- inch LCD

## Technical Specifications

Areas of Compression	Thigh, calf, and foot
Compression Modes	Uniform and sequential
Compression Pressure	Calf/Thigh: 20, 30, 40, 50, 60mmHg Foot: 120, 130, 140mmHg
Interval time	24, 48 and 60 seconds (selectable)

## Safety System

Alarms	Low Pressure, High Pressure, Sleeve not Connected, Battery Low, Battery Depleted, System Error
--------	--

Ingress Protection	IPX2
--------------------	------

## Power Supply

Mains Power	100-240 V AC, 50-60 Hz
Battery type	Li-ion battery (rechargeable)
Battery Life	Upto 8 hrs

## Accessories

	Power cord, Extension tube Sleeves - Uniform (calf/thigh) Sequential calf/thigh, Foot
--	---

## Warranty

	1 year
--	--------

# De Thromb

**YOUR PARTNER IN VASCULAR HEALTH**



## Why De Thromb ?

Deep Vein Thrombosis (DVT) is a serious condition where blood clots form in deep veins, typically in the legs. If left untreated, DVT can lead to potentially life threatening complications such as Pulmonary Embolism (PE), Venous Thrombo Embolism (VTE) De Thromb is a cutting-edge solution designed to effectively address deep vein thrombosis (DVT)

### vTitan Corporation Pvt. Ltd.

New Woodlands Building, Estancia IT Park,  
Plot no :140,151 Vallancheri Village,  
Chengalpattu District,  
Tamil Nadu - 603202, India  
Toll Free Number: 18005721664

[www.vtitan.com](http://www.vtitan.com)  
[sales@vtitan.com](mailto:sales@vtitan.com)

© 2025 vTitan Corporation Pvt. Ltd.  
980-1105 Rev 03

### Dynamic Intermittent Pneumatic Compression



Discover our groundbreaking dynamic compression technology designed to prevent Venous Thromboembolism (VTE).

This feature guarantees accurate and efficient compression therapy, enhancing blood circulation and lowering the chances of blood clot formation.

### Customizable Therapy Settings



Tailored treatment to each patient's unique needs with customizable therapy settings. Selectable pressure levels and deflation times allow for therapy treatment plans, optimizing comfort and efficacy.

### Selectable Modes



Uniform and Sequential modes that help to suit the specific needs of your therapy.

### Track Operating Hours



Provision to log pump usage for a particular patient and reset it after therapy



**Audio &  
Visual alarms**



**Attachable &  
Detachable hooks**

*Our cutting-edge DVT Compression Sleeves are designed to provide you with the ultimate comfort, whether you're recovering from surgery, managing a medical condition, or seeking proactive prevention measures.*



### Comfortable and Breathable Fabric

Crafted from premium-quality, breathable materials, our compression sleeves ensure maximum comfort even during extended wear. Say goodbye to discomfort and hello to all-day confidence!



### Customized Fit for Everybody

Available in a range of sizes, our sleeves offer a customized fit for individuals of all shapes and sizes. Whether you need a snug fit for optimal compression or a more relaxed feel for everyday wear, we've got you covered.



### Versatile Design

Our compression sleeves are versatile and perfect for wearing under clothing. Whether you're working, traveling, or relaxing, they offer the support you need wherever you go.



### Easy Care and Maintenance

Our sleeves offer hassle-free care and maintenance, maximizing longevity for consistent comfort and effectiveness during use.